

Wheelchair Restraints & Docking Station Rating Guide (v2)

Introduction:

To follow is a guide showing a range of wheelchair restraints and docking stations which are commonly used on vehicles in New Zealand. The table lists the maximum weight of wheelchair that the restraints are approved for. Restraints that incorporate attachment points for a lap belt are rated for the additional load, and the wheelchair maximum rating does not reduce when a lap belt is used. Wheelchairs that have an approved integral lap belt must use appropriate rated restraints.

Wheelchair restraints:

AMF-Bruns www.AMF-Bruns-mobility.com/products/wheelchair-restraints/

Protektor
 Wheelchairs up to:
85 kg

Black or Silver Series
 Hook, Karibiner or buckle attachment



Protektor HD160
 Wheelchairs up to:
160 kg

Platinum Series
 Hook, Karibiner or buckle attachment











Handi-Care www.handicare.com

Handi-Secure
 Wheelchairs up to:
100 kg

Hook or buckle attachment
 Product discontinued 2017



Koller www.Koller.co.uk	
Many different options available. For ratings refer to Koller website.	
Q'Straint www.QStraint.com	
QRT – Standard QRT – Deluxe QRT – Max Wheelchairs up to: 85 kg	 <p>Separate bracket for the lap belt attachment pin</p>  <p>One or two adjustment knobs</p>
QER-4000 Electric reel front only Wheelchairs up to: 85 kg	
QRT-3 QRT-350e 150 kg	 <p>Integrated lap belt attachment pin</p> <p>QRT-3 has a single loop hook and single webbing compared to the QRT-5</p>
QRT-5 QRT-550e Wheelchairs up to: 200 kg	 <p>QRT-5 has a double loop hook and double webbing</p>
QRT360 Wheelchairs up to: 85 kg* *Tested to WC18 specification with lap belt attached to wheelchair	 

Schnierle <https://www.schnierle.de/wheelchair-restraint-systems.html>

Quick-Secure/Safe-Trans
 Wheelchairs up to:
85 kg



Electro reel front only
 Wheelchairs up to:
85 kg



Safe-Trans 'Heavy Load'
 Wheelchairs up to:
160 kg
 Note the heavy-duty karabiner and black intermediate plate



Sure-Lock www.Sure-Lok.com

Titan 800
 Wheelchairs up to:
85 kg






Titan AL700
 Wheelchairs up to:
85 kg



Titan FF600
 Wheelchairs up to:
85 kg



Lindberg	
<p>Standard Wheelchairs up to: 85 kg</p> <p>Discontinued</p>	
Unwin (now part of Braunability) www.Braunability.eu	
<p>Many different options from 85 kg to 200 kg.</p> <p>For ratings refer to website</p>	

Docking stations	
Unwin (now part of Braunability) www.Braunability.eu	
<p>Dock'n'Lock Wheelchairs up to: 200 kg</p>	
EZ-Lock www.ezlock.net/	
<p>EZ-Lock Wheelchairs up to: 173 kg</p> <p>Approved for any wheelchair listed on the official website</p>	
Dahl www.dahlengineering.dk/en/products/dahl-docking-systems/dahl-docking-station-mk-ii/	
<p>MK11 Wheelchairs up to: 150 kg with lap belt on wheelchair 200 kg with lap belt to vehicle</p>	

Q'Straint www.QStraint.com	
QLK-100 Wheelchairs up to: 85 kg	
QLK-110 Wheelchairs up to: 150 kg	
QLK-150 Wheelchairs up to: 200 kg	
Uni-Loc www.freemob.com.au	
Wheelchairs up to: 140 kg	
Auto Transform www.autotransform.co.nz/solutions/independent-driving/	
Autodock Wheelchairs up to: 157 kg	

Finally:

If you have any queries, require any further clarification or have information suitable for inclusion in this Information Sheet, contact the technical team [by email](#).

A printed or saved copy of this information sheet may become out of date – the [LVVTA website](#) holds the latest document version.